



**Disconnect:
5 feet above
the floor**

CPI Generator

**Recommended Room Size:
9ft wide x 8ft deep x 8ft high**

**Minimum room size requirement:
7.7ft wide x 6ft deep x 7.2ft high**



PROPRIETARY AND CONFIDENTIAL
THE INFORMATION CONTAINED IN THIS
DRAWING IS THE SOLE PROPERTY OF
MEDICATECH USA. ANY
REPRODUCTION IN PART OR AS A WHOLE
WITHOUT THE WRITTEN PERMISSION OF
MEDICATECH USA IS
PROHIBITED.

TITLE: CRV2000- Min. room requirement		
SIZE: D	DWG. NO.	REV
SCALE: 1:20 WEIGHT:		SHEET 1 OF 1