



MEDICATECH USA®

MedicatechUSA presents
EXPRIMER 17X17 Wireless DR Panel


- ✓ Increased Dose Efficiency and Image Quality
- ✓ Improved Patient Throughput and Fast Workflow
- ✓ High DQE / MTF Value
- ✓ TRUVIEW®ART, Advanced Reverse Filtering Technology
- ✓ User Friendly & Ergonomic Design



Model DRTECH-1717W



Reliable
2.5GHz
5GHz
Wifi Performance



Errorless
Wireless



WPCS



OLED Display




IP54



USB C type





Excellent Dose Reduction




Exceptional Diagnostic Clarity



Focus Scintillator Technology

FOCUS (Fluorescent Optical Csi Upgraded Structure) has more focused optical signal compared to conventional Csi. As a result, ultra low noise is achieved by FOCUS technology.



Excellent Image Capturing capacity for all radiographic procedures and packed in a well-designed and easy to use casing for many years of use. Ideal for any radiographic needs.

Fast & Efficient Workflow


EXPRIMER Advanced series allows fast and efficient workflow brought by its fast exam time, speeding workflow and decreasing repeats.

First Image: Min. 4 sec.



Rugged Design with Portability

EXPRIMER Advanced series' LAN cable with fixed protector increases the reliability and durability. Moreover, attachable handle allows carrying detector more conveniently.



Model	DRTECH-1717W
Csl / Gadox	Csl / Gadox
Active area	430 x 430 (mm) / 17 x 17 (inch)
Resolution	3,072 X 3,072
Pixel pitch	140 μ m
Dimension	460 X 460 X 15 mm (Protruding Area: 24.2 x 125.5mm)
Weight	4.3 kg
X-ray Sync	Lossless AED/AWC, Sync Trigger
Data Interface	Gigabit Ethernet (LAN Cable)
ADC	16-bit



VOYANCE®

By Mediatech USA®

An industry leading complete software suite that delivers a secure and cost-effective way to acquire, view and share DICOM Images. Voyance® is one of the easiest, most user friendly interfaces. Its modern design & attention to the user experience along with its and features & functions make it a great tool to every user.

Instant installation and auto-configuration

Easy, user-friendly Modern user interface

Interface consistency across acquisition and viewing platforms

Optimum Workflow designed for efficiency

ClairVoyance® Image Processing Algorithm

Supports a wide range of detectors and generators

X-RAY Generator integration (CPI & DRGEM)

Auto Stitching and Auto Cropping

Side by side comparisons

Fully DICOM 3.0 Compatible: MWL, MPPS, and all DICOM ops.

Reject Analysis





MEDICATECH USA®



T. +1 800-817-5030

T. +1 949-679-2881

F. +1 949-679-2882



www.medicattechusa.com



info@medicattechusa.com