



Recommended Room Size:
9ft wide x 8ft deep x 8ft high

Minimum room size requirement:

TITLE:
CRV1500- Min. room
requirement

7.7ft wide x 6ft deep x 7.2ft high

TOLERANCES UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN MILLIMETERS
TOLERANCES:
.XXX ±0.010 ANGLE: ±0.5'
.XX ±0.05
.X ±0.2
SURFACE FINISH: 0.8-1.6µm

Projection Third Angle

DRAWN USMedimade	1/11/2022
APPROVED	
This drawing and any information or descriptive materials set out on it, are the confidential and copyright property of USMedimade, INC and must not be disclosed, copied, loaned in whole or part or used for any purpose without the written permission of USMedimade, INC.	
PART #	CR2 TABLE ASSY

TITLE CRV 1500			
SIZE B	Material	DWG NO CRV 1500.idw	REV
SCALE 1/17	Finish	SHEET 1 OF 1	