

DOG

Body Part	Size	kVp	mA	s	mAs
SKULL / SKULL LAT	1-4 cm	72	320	0.01	3.2
SKULL / SKULL LAT	5-8 cm	77	320	1.20E-02	4
SKULL / SKULL LAT	9-12 cm	73	320	1.60E-02	5
SKULL / SKULL LAT	13-16 cm	77	320	0.02	6.4
SKULL / SKULL LAT	17-20 cm	82	320	2.50E-02	8
SKULL / SKULL LAT	21-24 cm	75	320	2.50E-02	8
SKULL / SKULL LAT	25-30 cm	75	320	2.50E-02	8

SKULL / SKULL VD	1-4 cm	72	320	0.01	3.2
SKULL / SKULL VD	5-8 cm	77	320	1.20E-02	4
SKULL / SKULL VD	9-12 cm	77	320	1.60E-02	5
SKULL / SKULL VD	13-16 cm	77	320	0.02	6.4
SKULL / SKULL VD	17-20 cm	77	320	2.50E-02	8
SKULL / SKULL VD	21-24 cm	80	320	2.50E-02	8
SKULL / SKULL VD	25-30 cm	80	320	2.50E-02	8

SKULL / SKULL DV	1-4 cm	72	320	0.01	3.2
SKULL / SKULL DV	5-8 cm	77	320	1.60E-02	5
SKULL / SKULL DV	9-12 cm	72	320	2.50E-02	8
SKULL / SKULL DV	13-16 cm	77	320	0.02	6.4
SKULL / SKULL DV	17-20 cm	77	320	3.20E-02	10
SKULL / SKULL DV	21-24 cm	80	320	0.04	12.5
SKULL / SKULL DV	25-30 cm	80	320	2.50E-02	8

R.EXT / R.EXT LAT	1-4 cm	72	250	0.01	2.5
R.EXT / R.EXT LAT	5-8 cm	72	320	0.01	3.2
R.EXT / R.EXT LAT	9-12 cm	72	320	1.20E-02	4
R.EXT / R.EXT LAT	13-16 cm	77	320	1.60E-02	5
R.EXT / R.EXT LAT	17-20 cm	72	320	0.02	6.4
R.EXT / R.EXT LAT	21-24 cm	75	320	2.50E-02	8
R.EXT / R.EXT LAT	25-30 cm	75	320	2.50E-02	8

R.EXT / R.EXT VD	1-4 cm	72	250	0.01	2.5
R.EXT / R.EXT VD	5-8 cm	72	320	0.01	3.2
R.EXT / R.EXT VD	9-12 cm	72	320	1.20E-02	4

R.EXT / R.EXT VD	13-16 cm	77	320	1.60E-02	5
R.EXT / R.EXT VD	17-20 cm	72	320	0.02	6.4
R.EXT / R.EXT VD	21-24 cm	77	320	2.50E-02	8
R.EXT / R.EXT VD	25-30 cm	80	320	0.04	12.5

R.EXT / R.EXT DV	1-4 cm	72	250	0.01	2.5
R.EXT / R.EXT DV	5-8 cm	72	320	0.01	3.2
R.EXT / R.EXT DV	9-12 cm	72	320	1.20E-02	4
R.EXT / R.EXT DV	13-16 cm	77	320	1.60E-02	5
R.EXT / R.EXT DV	17-20 cm	72	320	0.02	6.4
R.EXT / R.EXT DV	21-24 cm	77	320	2.50E-02	8
R.EXT / R.EXT DV	25-30 cm	80	320	2.50E-02	8

F.EXT / F.EXT LAT	1-4 cm	72	250	0.01	2.5
F.EXT / F.EXT LAT	5-8 cm	72	320	0.01	3.2
F.EXT / F.EXT LAT	9-12 cm	72	320	1.20E-02	4
F.EXT / F.EXT LAT	13-16 cm	77	320	1.60E-02	5
F.EXT / F.EXT LAT	17-20 cm	82	320	0.02	6.4
F.EXT / F.EXT LAT	21-24 cm	85	320	2.50E-02	8
F.EXT / F.EXT LAT	25-30 cm	85	320	2.50E-02	8

F.EXT / F.EXT VD	1-4 cm	72	250	0.01	2.5
F.EXT / F.EXT VD	5-8 cm	72	320	0.01	3.2
F.EXT / F.EXT VD	9-12 cm	72	320	1.20E-02	4
F.EXT / F.EXT VD	13-16 cm	77	320	1.60E-02	5
F.EXT / F.EXT VD	17-20 cm	82	320	0.02	6.4
F.EXT / F.EXT VD	21-24 cm	80	320	3.20E-02	10
F.EXT / F.EXT VD	25-30 cm	80	320	2.50E-02	8

F.EXT / F.EXT DV	1-4 cm	72	250	0.01	2.5
F.EXT / F.EXT DV	5-8 cm	72	320	0.01	3.2
F.EXT / F.EXT DV	9-12 cm	72	320	1.20E-02	4
F.EXT / F.EXT DV	13-16 cm	77	320	1.60E-02	5
F.EXT / F.EXT DV	17-20 cm	72	320	0.02	6.4
F.EXT / F.EXT DV	21-24 cm	77	320	0.02	6.4
F.EXT / F.EXT DV	25-30 cm	80	320	2.50E-02	8

THORAX / THORAX LAT	1-4 cm	72	160	0.01	1.6
THORAX / THORAX LAT	5-8 cm	77	200	0.01	2
THORAX / THORAX LAT	9-12 cm	82	200	0.01	2
THORAX / THORAX LAT	13-16 cm	87	200	0.01	2
THORAX / THORAX LAT	17-20 cm	100	200	0.01	2
THORAX / THORAX LAT	21-24 cm	112	200	0.01	2
THORAX / THORAX LAT	25-30 cm	110	250	0.01	2.5

THORAX / THORAX VD	1-4 cm	72	160	0.01	1.6
THORAX / THORAX VD	5-8 cm	77	200	0.01	2
THORAX / THORAX VD	9-12 cm	82	200	0.01	2
THORAX / THORAX VD	13-16 cm	90	200	0.01	2
THORAX / THORAX VD	17-20 cm	100	200	0.01	2
THORAX / THORAX VD	21-24 cm	112	200	0.01	2
THORAX / THORAX VD	25-30 cm	110	250	0.01	2.5

THORAX / THORAX DV	1-4 cm	72	160	0.01	1.6
THORAX / THORAX DV	5-8 cm	77	200	0.01	2
THORAX / THORAX DV	9-12 cm	82	200	0.01	2
THORAX / THORAX DV	13-16 cm	90	200	0.01	2
THORAX / THORAX DV	17-20 cm	100	200	0.01	2
THORAX / THORAX DV	21-24 cm	112	200	0.01	2
THORAX / THORAX DV	25-30 cm	110	250	1.20E-02	3.2

SPINE / SPINE LAT	1-4 cm	72	320	0.01	3.2
SPINE / SPINE LAT	5-8 cm	72	320	1.20E-02	4
SPINE / SPINE LAT	9-12 cm	72	320	1.60E-02	5
SPINE / SPINE LAT	13-16 cm	77	320	0.02	6.4
SPINE / SPINE LAT	17-20 cm	77	320	2.50E-02	8
SPINE / SPINE LAT	21-24 cm	77	320	3.20E-02	10
SPINE / SPINE LAT	25-30 cm	82	320	0.04	12.5

SPINE / SPINE VD	1-4 cm	72	320	0.01	3.2
SPINE / SPINE VD	5-8 cm	72	320	1.20E-02	4
SPINE / SPINE VD	9-12 cm	72	320	1.60E-02	5

SPINE / SPINE VD	13-16 cm	72	320	0.02	6.4
SPINE / SPINE VD	17-20 cm	77	320	2.50E-02	8
SPINE / SPINE VD	21-24 cm	82	320	3.20E-02	10
SPINE / SPINE VD	25-30 cm	82	320	0.05	16

SPINE / SPINE DV	1-4 cm	66	320	3.20E-02	10
SPINE / SPINE DV	5-8 cm	69	320	3.20E-02	10
SPINE / SPINE DV	9-12 cm	70	320	0.04	12.5
SPINE / SPINE DV	13-16 cm	75	320	0.04	12.5
SPINE / SPINE DV	17-20 cm	80	320	0.04	12.5
SPINE / SPINE DV	21-24 cm	77	320	0.05	16
SPINE / SPINE DV	25-30 cm	82	320	0.05	16

ABDOMEN / ABDOM LT	1-4 cm	72	320	0.01	3.2
ABDOMEN / ABDOM LT	5-8 cm	77	320	1.20E-02	4
ABDOMEN / ABDOM LT	9-12 cm	82	320	1.20E-02	4
ABDOMEN / ABDOM LT	13-16 cm	77	320	1.60E-02	5
ABDOMEN / ABDOM LT	17-20 cm	77	320	0.02	6.4
ABDOMEN / ABDOM LT	21-24 cm	82	320	0.02	6.4
ABDOMEN / ABDOM LT	25-30 cm	82	320	2.50E-02	8

ABDOMEN / ABDOM VD	1-4 cm	72	320	0.01	3.2
ABDOMEN / ABDOM VD	5-8 cm	77	320	1.60E-02	5
ABDOMEN / ABDOM VD	9-12 cm	77	320	1.20E-02	4
ABDOMEN / ABDOM VD	13-16 cm	77	320	1.60E-02	5
ABDOMEN / ABDOM VD	17-20 cm	77	320	0.02	6.4
ABDOMEN / ABDOM VD	21-24 cm	82	320	2.50E-02	8
ABDOMEN / ABDOM VD	25-30 cm	82	320	2.50E-02	8

ABDOMEN / ABDOM DV	1-4 cm	72	320	1.60E-02	5
ABDOMEN / ABDOM DV	5-8 cm	77	320	1.60E-02	5
ABDOMEN / ABDOM DV	9-12 cm	77	320	0.02	6.4
ABDOMEN / ABDOM DV	13-16 cm	77	320	2.50E-02	8
ABDOMEN / ABDOM DV	17-20 cm	82	320	2.50E-02	8
ABDOMEN / ABDOM DV	21-24 cm	82	320	3.20E-02	10
ABDOMEN / ABDOM DV	25-30 cm	87	320	3.20E-02	10

PELVIS / PELVIS LAT	1-4 cm	72	320	0.01	3.2
PELVIS / PELVIS LAT	5-8 cm	72	320	1.20E-02	4
PELVIS / PELVIS LAT	9-12 cm	77	320	1.60E-02	5
PELVIS / PELVIS LAT	13-16 cm	77	320	0.02	6.4
PELVIS / PELVIS LAT	17-20 cm	82	320	0.02	6.4
PELVIS / PELVIS LAT	21-24 cm	82	320	2.50E-02	8
PELVIS / PELVIS LAT	25-30 cm	82	320	3.20E-02	10

PELVIS / PELVIS VD	1-4 cm	72	320	0.01	3.2
PELVIS / PELVIS VD	5-8 cm	72	320	1.20E-02	4
PELVIS / PELVIS VD	9-12 cm	77	320	1.60E-02	5
PELVIS / PELVIS VD	13-16 cm	77	320	0.02	6.4
PELVIS / PELVIS VD	17-20 cm	82	320	0.02	6.4
PELVIS / PELVIS VD	21-24 cm	82	320	2.50E-02	8
PELVIS / PELVIS VD	25-30 cm	82	320	3.20E-02	10

PELVIS / PELVIS DV	1-4 cm	72	320	2.50E-02	8
PELVIS / PELVIS DV	5-8 cm	77	320	2.50E-02	8
PELVIS / PELVIS DV	9-12 cm	77	320	1.20E-02	4
PELVIS / PELVIS DV	13-16 cm	77	320	1.60E-02	5
PELVIS / PELVIS DV	17-20 cm	82	320	0.02	6.4
PELVIS / PELVIS DV	21-24 cm	82	320	2.50E-02	8
PELVIS / PELVIS DV	25-30 cm	87	320	0.04	12.5